

CA190

TILTMA

Tilt Table (Electric)



- 1 Year on Gas Struts, Castors, Upholstery
- 2 Years on Electronic Components
- 3 Years on Frame & Pillar



Prussian Blue Licorice Black Dim Grey Merlot (Red Wine) Creamy Beige



2 Motors



Hand Control

Maximum User Weight
225 kgs / 496 lbs

Product Features

This is a CA190 Tiltma and it is the Ultima of Tilt Tables for diagnosis of syncope and/or dysautonomia. In terms of its comfort rating, mechanical stability and operational movement, it is one of the best. The two adjustable patient's foot platform (Foot Paddle) can be tilted to a positive angle of 20° degree.

The CA190 Tiltma is equipped with the following accessories in the package:

- Chest Strap * 1
- Thigh Straps * 2
- Calf Straps * 2
- Working Table * 1

All 4 corners of the CA190 Tiltma are supported by a 5" inch castor with Swivel Lock and Brake. The Frame Construction is made with Steel Tubes of (a) 6 * 3cm and (b) 7 * 3cm Dimensions and by the thickness of 2mm.

For absolute safety requirement, the CA190 Tiltma is equipped with a "1-PULL EMERGENCY HANDLE" at the back of table. When in an emergency situation where no utility power is available, the emergency handle can be pulled to allow the table top to be moved back from vertical to horizontal position and the patient can then be nstrapped and transferred to normal patient's bed.

In any given case of loss of utility power, just pull on the RED EMERGENCY HANDLE at the back of table to disengage the Motorized Actuator (otherwise locked) to allow the patient to be unstrapped and lifted off the table.

CA190 Tiltma is equipped with 2 motors and controlled by a (2 Buttons) Hand Control:

- 1st Motorized Actuator: Height Adjustment (55~90cm)
- 2nd Motorized Actuator: Tilt Angle Adjustment (0°~90°)


There is 1 set of UP/DOWN Buttons on the Hand Control to control both Elevation and Tilting.

Up Button: When the table elevates to the highest position(90cm/35.4in), keep pressing the button and the tilt function will be activated.

Down Button: When the table tilts back to the horizontal position, keep pressing the button and the table will be lowering down.



1-Pull Emergency Handle



Technical Specifications

