

DS MAREF
DAESUNG MAREF CO.,LTD.



CTC-7

Cryo.Thermo.Compression system



Specifications

- + Power source : AC 100-240V, 50/60Hz
- + Power consumption : 270VA
- + Mode : Cold, Heat, Contrast
- + Temperature : Cold : 6°C to 15°C ±3°C
Heat : 35°C to 45°C ±3°C
- + Time Range : Cold 5 - 60 min / Preset mode
Heat 5 - 30 min / Preset mode
- + Dimensions : 27.2(W) x 18(D) x 23.9(H) cm
- + Weight : 3.7kg / 8.2lb (Empty water container)
3.9kg / 8.6lb (Full water container)
- + Type of protection IP : IPX0
- + Recommended Liquid : Distilled water
- + Class / type of protection : Class I/BF-type
- + SAFETY : IEC 60601-1
- + EMC : IEC 60601-1-2

The CTC-7 is an all in one rehabilitation system allowing you to customise temperature, compression level and cycle time for optimum treatment.

Unlike Ice-refilling systems, that fluctuate in temperature as the ice melts, CTC-7 offers continuous treatment to the therapy area with a consistent temperature for maximum effectiveness.

Iceless & Continuous Therapy

Relieving Pain

Reducing Inflammation

Reducing Swelling



Faster Recovery Time

Cost Savings



The CTC-7 Wraps

are designed to maximize the functionality of the CTC-7 system by providing cold or hot water with air compression.

The compression pad covers the inner water pad so that cold/hot water can flow evenly over the treatment area, with the water pad being closer to the skin. Alongside periodic air compression to make cold/heat treatment more effective. In addition, customised compression levels and interval times can be applied intermittently to best suit treatment needs.

Appropriate straps are attached to the wrap so that it can be adjusted and secured to fit a variety of sizes.

CTC wrap has its own drain valve, so it is easy to keep it dry after draining it when storing it after use.



Key functions

Features

Iceless
Operation

Continuous
Use Time

Superior
Wrap Design

Protocols

Cold
+ Compression

Heat
+ Compression

Contrast
+ Compression

Convenience

Auto
Water Check &
Alarm

Steam
Cleaning

Auto
Drain

Single Person Use Wraps



CTCAFD0



CTCUMD0 · CTCULD0



CTCKFD0



CTCEF00



CTCSMD0



CTCFE00

Order No.	One Person Use Wraps List	Size
CTCKFD0	KNEE WRAP	UNIVERSAL
CTCUMD0	UNIVERSAL WRAP	M
CTCULD0	UNIVERSAL WRAP	L
CTCBFD0	BACK WRAP	UNIVERSAL
CTCEF00	ELBOW WRAP	UNIVERSAL
CTCSMD0	SHOULDER WRAP	M
CTCAFD0	ANKLE WRAP	UNIVERSAL

Re-usable Wraps



CTCAF00



CTCUM00 · CTCUM00



CTCKF00



CTCEF00



CTCSM00



CTCFE00

Order No.	Resuable Wraps List	Size
CTCKF00	KNEE WRAP	UNIVERSAL
CTCUM00	UNIVERSAL WRAP	M
CTCUL00	UNIVERSAL WRAP	L
CTCBF00	BACK WRAP	UNIVERSAL
CTCEF00	ELBOW WRAP	UNIVERSAL
CTCSM00	SHOULDER WRAP	M
CTCAF00	ANKLE WRAP	UNIVERSAL

DS MAREF
DAESUNG MAREF CO.,LTD.

Maydwell Avenue, Stane Street, Slinfold,
Horsham, West Sussex, RH13 0GN

Tel: 01403 799190

E-mail: medsales@gandn.com

www.gandn.com

