

Physiotherapy & Rehabilitation



YouTecarX



CHOOSE THE IDEAL CONFIGURATION FOR YOUR NEEDS AND THOSE OF YOUR PATIENTS

WHAT IS TECAR THERAPY? WHAT BENEFIT DOES IT OFFER?

CONTACT DIATHERMY – more commonly known as TECAR THERAPY – is a medical practice that applies radiofrequency to biological tissue with a frequency ranging from 300 kHz to 1 MHz without generating muscle contraction.

The position of the electrodes on the body determines the current path, thereby ensuring a deep, targeted, comfortable application that is also easy to control and manage.

The heat-induced vasodilation increases the exchange of substances promoting the drainage of the inflamed area and the influx of reparative cells.

Tecartherapy is based on the use of an electromedical device that uses high-frequency currents for the **CHEMICAL | THERMAL | MECHANICAL** biostimulation of tissues with different resistivity.

Physiotherapy, rehabilitation, and pain management are the fields that most use the beneficial effects of tecartherapy.

Tecar Therapy units are often the first piece of equipment to be bought.

THE NUMBER OF UNITS IS GROWING ALL OVER THE WORLD!

WHAT FIELD DO YOU WORK IN?

- Sports Rehabilitation
- Dermal-Functional Rehabilitation
- Physio-Aesthetics
- Pelvic Floor Rehabilitation
- MuskuloSkeletal Rehabilitation:
 - Cronic & Acute
 - Post-Traumatic
 - Post-Surgical

WHAT TECHNIQUES DO YOU COMBINE?

- Fascia Therapy
- Manual Therapy
- Therapeutic Exercise
- Osteopathy



CONFIGURATION STEP BY STEP

CREATE YOUR OWN TECAR THERAPY DEVICE

BASE CONFIGURATION:


Power: 150 Watt
Frequency: 500 KHz
Output: Unipolar
Emission: Continuous

INCLUDES:

- 1 power cable
- 1 cable for unipolar applications
- 1 unipolar handpiece
- 1 steel neutral electrode
- 1 resistive electrode Ø 61 mm
- 1 capacitive electrode Ø 61 mm
- Cable with clamp for disposable neutral electrode
- Disposable gel adhesive electrode (5 pcs)



CUSTOMIZATION OPTIONS



THE PRICE? YOU CHOOSE!



BASE CONFIGURATION

CUSTOMIZATIONS

TOP MODEL

€

€€

€€€

€€€€€

YOULECAR



PAGE 4

CORE & SYSTEM

Choose the power, frequencies, and forms of emission



CHOOSE THE POWER

- 150 WATT (included)
- 200 WATT
- 250 WATT




CHOOSE THE FREQUENCIES

- 500 KHZ (included)
- STANDARD PACK 3 FREQ. (500 KHz / 750 KHz / 1 MHz)
- ON DEMAND PACK (500 KHz + a choice of three frequencies from 300 kHz to 1 MHz)



CHOOSE THE EMISSION

- CONTINUOUS (included)
- SENSITIVE ATHERMY
- HEA - HIGH ENERGY ATHERMY
- ATHERMY PACK (SENSITIVE + HEA)




AUX OUTPUT

The doorway to a whole world of new treatment possibilities

BOOST THE POWER OF YOUTECARX!



The multifunction AUX output opens the door to the highest performance levels using special accessories and application systems.

JOULE MODE is an innovative operating method that includes a dedicated hardware and software system



ACCESSORIES

Applicators designed to reach specific therapeutic targets



UNIPOLAR

Complete your accessory kit for traditional unipolar mode operation
(active electrode + neutral plate)




AUTOMATIC APPS

Not operator-dependent




NEUTRO DINAMICO®

Dynamic applications with two active handpieces and without neutral plate, using NeutroDinamico® accessories



BIPOLAR

Work without neutral plate using
Bipolar accessories



INTIMITY

Pelvic Floor rehabilitation




TRIGGER & FASCIA



PHYSIO-AESTHETICS

Physio-aesthetics Accessories



OCTOPUS CONTOUR

Extreme freedom for every kind of treatment



YOU1ECAR

PAGE 15



cnergy[®]

**ONE STOP SHOP
FOR ALL YOUR
PHYSIO AND REHAB
NEEDS**

Bracing & Supports
Physiotherapy Modalities
rPMS

Spinal Decompression Therapy
Transcranial Magnetic Stimulation

Hot & Cold Therapy

Exercise Therapy

Cryo Chamber

Treatment Tables

Hydrotherapy Pool

Underwater Treadmill

Virtual Rehab

Balancing Platform

GAIT Lab

Isokinetic System

NeuroRehabilitation Robot

Robotic GAIT Trainer

Spine & Whole Body Analysis

Imported & Marketed by

Cnergy Incorporation

5277 | 8th Street | Ram Nagar North Extension | Madipakkam |
Chennai 600091 | India, T: +91 44 4006 7974, +91 72000 53696 |
incorporation@cnergy.in | customercare@cnergy.in | www.cnergy.in

Manufactured by
Marpe Co., Ltd
South Korea